



Native Rat Thrombin (DAG-WT3315)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Rat
Purity	> 90% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Format	Liquid
Size	100 µg
Buffer	50 mM sodium citrate, 0.2 M NaCl, pH 6.5
Preservative	0.1% PEG-8000
Storage	Store at -80°C

BACKGROUND

Keywords	Thrombin; alpha thrombin; Coagulation factor II
----------	-------------------------------------------------