



Native Bovine α Thrombin (DAG-WT3316)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Bovine
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	37 kDa
Format	Liquid
Size	1000 units
Buffer	50 mM sodium citrate, 0.2 M NaCl, pH 6.5
Preservative	0.1% PEG-8000
Storage	Store at -80°C

BACKGROUND

Keywords	Thrombin; alpha thrombin; Coagulation factor II
----------	---