



# Native Bovine Prothrombin (DAG-WT3317)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Bovine
Purity	> 90% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	70 kDa
Format	Liquid
Size	10 mg
Buffer	20 mM Tris-HCl, 0.1 M NaCl, pH 7.4
Storage	Store at -80°C

## BACKGROUND

Keywords	PIVKA-II; Prothrombin; DCP
----------	----------------------------