



Native Bovine Antithrombin (DAG-WT3318)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Bovine
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	58 kDa
Format	Liquid
Size	1 mg
Buffer	20 mM Tris-HCl, 0.1M Sodium Citrate, 0.15 M NaCl, pH 8.3
Storage	Store at -80°C

BACKGROUND

Keywords	Antithrombin; AT; ATIII; SERPINC1
----------	-----------------------------------