



Recombinant Timothy Grass Phl p 5a [His, Biotin] (DAG-WT332)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Timothy Grass Phl p 5 has been expressed as a recombinant antigen fused to a his-tag in its N-terminus and in vivo biotinylated
Purity	> 95% , as determined by SDS-PAGE
Conjugate	His, Biotin
Applications	WB/ELISA
Molecular Weight	40.3 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	20 mM phosphate buffer pH 6.0 and 0.15 M NaCl and 0.1% polyoxyethylene (10) tridecyl ether
Preservative	None
Storage	Store at -20°C to -80°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords Timothy Grass; Phl p 5; Allergen