



Native Canine Fibrinogen (DAG-WT3320)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Canine
Purity	> 90% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Format	Liquid
Size	1 mg
Buffer	20 mM sodium citrate-HCl, pH 7.4
Storage	Store at -80°C

BACKGROUND

Keywords Fibrinogen; Factor I; Fibrinogen alpha chain; Fibrinogen beta chain; Fibrinogen gamma chain