



# Native Murine Fibrinogen (DAG-WT3321)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Murine
Purity	> 90% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Format	Liquid
Size	1 mg
Buffer	20 mM sodium citrate-HCl, pH 7.4
Storage	Store at -80°C

## BACKGROUND

Keywords	Fibrinogen; Factor I; Fibrinogen alpha chain; Fibrinogen beta chain; Fibrinogen gamma chain
----------	---