



## Native Rat Fibrinogen (DAG-WT3323)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Species	Rat
Purity	> 90% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Format	Liquid
Size	1 mg
Buffer	20 mM sodium citrate-HCl, pH 7.4
Storage	Store at -80°C

### BACKGROUND

**Keywords** Fibrinogen; Factor I; Fibrinogen alpha chain; Fibrinogen beta chain; Fibrinogen gamma chain