



Native Human Plasmin (DAG-WT3325)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 90% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	83 kDa
Format	Liquid
Size	1 mg
Buffer	50 mM Hepes, 50 mM Sodium Acetate, pH 8.5, 50% Glycerol
Storage	Store at -80°C

BACKGROUND

Keywords	PLG; Plasminogen; Plasmin; HAE4
----------	---------------------------------