



# Native Human Protein S (DAG-WT3326)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	69 kDa
Format	Liquid
Size	1 mg
Buffer	20 mM Tris-HCl, 0.1 M NaCl, 1 mM Benzamidine, pH 7.4
Storage	Store at -80°C

## BACKGROUND

Keywords	Protein S; PROS; PROS1; Coagulation
----------	-------------------------------------