



## Native Human Protein Z (DAG-WT3327)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Format	Liquid
Size	500 µg
Buffer	50 mM Tris-HCl, 0.1 M NaCl, pH 7.5
Storage	Store at -80°C

### BACKGROUND

Keywords	Protein Z; PZ; PROZ; Coagulation
----------	----------------------------------