



Native Alpha2-Antiplasmin (DAG-WT3328)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	67 kDa
Format	Liquid
Size	100 µg
Buffer	20 mM Tris-HCl, 0.15 M NaCl, pH 7.4
Storage	Store at -80°C

BACKGROUND

Keywords	Alpha 2-Antiplasmin; α 2-Antiplasmin; α2AP; α2 Plasmin inhibitor
----------	--