



Native Corn Trypsin Inhibitor (DAG-WT3329)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 90% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Format	Liquid
Size	1 mg
Buffer	20 mM Tris-HCl, 0.3 M NaCl, pH 8.2
Storage	Store at -80°C

BACKGROUND

Keywords	Trypsin inhibitor; TI; Serine protease inhibitor; Serpin
----------	--