



# Native Human Glycoprotein IIb/IIIa (DAG-WT3330)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Human
Conjugate	N/A
Applications	Immunoassays
Format	Liquid
Size	1 mg
Buffer	20 mM Tris-HCl, 0.1 M NaCl, 1 mM CaCl <sub>2</sub> , 50% glycerol, pH 7.4
Preservative	0.05% NaN <sub>3</sub>
Storage	Store at -80°C

## BACKGROUND

Keywords	GP1Ib111a; Glycoprotein 11b/111a; GP11b; GP111a
----------	---