



Native PIIINP Antigen (DAG-WT3349)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
---------------	---------------------------------

Conjugate	N/A
------------------	-----

Applications	Immunoassays
---------------------	--------------

Format	Liquid
---------------	--------

Size	1 mg
-------------	------

Storage	Store at -80°C
----------------	----------------

BACKGROUND

Keywords	PIIINP; N-terminal propeptide of procollagen type III
-----------------	---