



## Native AMH Antigen (DAG-WT3353)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|              |                                 |
|--------------|---------------------------------|
| Purity       | > 95% as determined by SDS-PAGE |
| Conjugate    | N/A                             |
| Applications | Immunoassays                    |
| Format       | Liquid                          |
| Size         | 1 mg                            |
| Storage      | Store at -80°C                  |

### BACKGROUND

|          |                             |
|----------|-----------------------------|
| Keywords | AMH; Anti-Müllerian Hormone |
|----------|-----------------------------|