



# Recombinant DENV type 3 Envelope Protein (DAG-WT3360)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	The DENV E protein consists of three domains (designated E protein domain EDI, EDII and EDIII)
Purity	> 90% as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	ELISA, CLIA, LFIA
Molecular Weight	43 kDa
Format	Liquid
Size	100 µg, 1 mg
Buffer	10 mM PBS, pH7.4
Storage	Store at -20°C

## BACKGROUND

Keywords	Dengue virus; DENV; Envelope protein; E protein
----------	---