



# Recombinant DENV type 1 NS2 Protein (DAG-WT3363)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	ELISA, CLIA, LFIA
Format	Liquid
Size	100 µg, 1 mg
Buffer	10 mM PBS, pH7.4
Storage	Store at -20°C

## BACKGROUND

Keywords	Dengue virus; DENV; Nonstructural protein 1
----------	---