



# Recombinant Human PIVKA-II Protein [His] (DAG-WT3369)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

**Purity** > 90% as determined by SDS-PAGE

**Conjugate** His

**Applications** ELISA, CLIA

**Format** Liquid

**Size** 1 mg

**Buffer** 20 mM Tris, 150 mM NaCl, pH8.0

**Storage** Store at -20°C

## BACKGROUND

**Keywords** PIVKA-II; Prothrombin; DCP