



Recombinant CMV gB (UL55) (aa 549-645) [GST] (DAG-WT337)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% , as determined by SDS-PAGE
Conjugate	GST
Applications	WB/ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	100 µg
Buffer	1.5 M Urea; 25 mM Tris HCl pH 8.0; 0.2% Triton X, 50% glycerol
Preservative	None
Storage	Store at -80°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Introduction	Cytomegalovirus (from the Greek cyto-, "cell", and -megal-, "large") is a viral genus of the viral family known as Herpesviridae or herpesviruses. It is typically abbreviated as CMV. The species that infects humans is commonly known as human CMV (HCMV) or human herpesvirus-5 (HHV-5), and is the most studied of all cytomegaloviruses.
Keywords	CMV; CMV glycoprotein B; gB; Glycoprotein B; Glycoprotein GP55; HCMV; HHV5 glycoprotein B; UL55; CMV gB