



Recombinant *Mytilus edulis* Myt e 1 [His] (DAG-WT3372)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA, CLIA, WB
Format	Liquid
Size	1 mg
Buffer	20 mM Tris, 150 mM NaCl, pH8.0
Storage	Store at -20°C

BACKGROUND

Keywords	<i>Mytilus edulis</i> ; Myt e 1
----------	---------------------------------