



# Recombinant *Blattella germanica* Bla g 7 (DAG-WT3435)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	33.8 kDa
Format	Lyophilized
Size	250 µg
Buffer	Lyophilized from sterile 20 mM Tris, 300 mM NaCl, 200 mM Arg, pH 7.4
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	<i>Blattella germanica</i> ; Bla g 7
----------	--------------------------------------