



Recombinant *Canis familiaris* Can f 4 (DAG-WT3437)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	19.1 kDa
Format	Lyophilized
Size	250 µg
Buffer	Lyophilized from sterile 50 mM Tris, 100 mM NaCl, 0.02 % TWEEN 80, pH 8.0
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	<i>Canis familiaris</i> ; Can f 4
----------	-----------------------------------