



Recombinant Canis familiaris Can f 6 (DAG-WT3438)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	21.6 kDa
Format	Lyophilized
Size	250 µg
Buffer	Lyophilized from sterile PBS, pH 7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Canis familiaris; Can f 6
----------	---------------------------