



# Recombinant *Aspergillus fumigatus* Asp f 22 (DAG-WT3439)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	48.1 kDa
Format	Lyophilized
Size	250 µg
Buffer	Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 30% Glycerol, pH 7.5
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	<i>Aspergillus fumigatus</i> ; Asp f 22
----------	---