



Recombinant DENV Nonstructural Protein 1n [His] (DAG-WT344)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% , as determined by SDS-PAGE
Conjugate	His
Applications	WB/ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	100 µg
Buffer	100 mM Tris HCl pH 8.0; 8M Urea, 200mM NaCl
Preservative	None
Storage	Store at -80°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Introduction	NS1 is one of 7 Dengue Virus non-structural proteins which are thought to be involved in viral replication. NS1 exists as a monomer in its immature form but is rapidly processed in the endoplasmic reticulum to form a stable dimer. A small amount of NS1 remains associated with intracellular organelles where it is thought to be involved in viral replication. The rest of NS1 is found either associated with the plasma membrane or secreted as a soluble hexadimer. NS1 is essential for viral viability but its precise biological function is unknown. Antibodies raised in response to NS1 in viral infection can cross react with cell surface antigens on epithelial cells and platelets and this has been implicated in the development of Dengue Hemorrhagic fever.
---------------------	---

Keywords

Dengue Virus; Non-structural protein; NS1
