



Recombinant *Aspergillus fumigatus* Asp f 28 (DAG-WT3441)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	12.9 kDa
Format	Lyophilized
Size	250 µg
Buffer	Lyophilized from sterile 20 mM Tris, 300 mM NaCl, pH 7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	<i>Aspergillus fumigatus</i> ; Asp f 28
----------	---