



Recombinant *Cupressus arizonica* Cup a 1 (DAG-WT3445)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	38.6 kDa
Format	Lyophilized
Size	250 µg
Buffer	Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 200 mM Arg, 2 M Urea, pH 8.0
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Cupressus arizonica; Cup a 1
----------	------------------------------