



Recombinant *Ambrosia artemisiifolia* Amb a 11 (DAG-WT3446)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	41.7 kDa
Format	Lyophilized
Size	250 µg
Buffer	Lyophilized from sterile 20 mM Tris, 300 mM NaCl, 200 mM Arg, 2 M Urea, pH 7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	<i>Ambrosia artemisiifolia</i> ; Amb a 11
----------	---