



# Recombinant *Amaranthus retroflexus* Ama r 2 (DAG-WT3447)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	15.2 kDa
Format	Lyophilized
Size	250 µg
Buffer	Lyophilized from sterile 20 mM Tris, 200 mM NaCl, pH 7.4
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	<i>Amaranthus retroflexus</i> ; Ama r 2
----------	---