



# Recombinant *Parietaria judaica* Par j 2 (DAG-WT3448)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	12.3 kDa
Format	Lyophilized
Size	250 µg
Buffer	Lyophilized from sterile PBS, pH 7.4
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	<i>Parietaria judaica</i> ; Par j 2
----------	-------------------------------------

## GENE INFORMATION

UniProt ID	P55958
------------	--------