



# Recombinant *Arachis hypogaea* Ara h 1 (DAG-WT3453)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	70.1 kDa
Format	Lyophilized
Size	250 µg
Buffer	Lyophilized from sterile 20 mM Tris, 300 mM NaCl, 10% Glycerol, pH 7
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	<i>Arachis hypogaea</i> ; Ara h 1
----------	-----------------------------------