



# Recombinant *Bertholletia excelsa* Ber e 1 (DAG-WT3456)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	14.2 kDa
Format	Lyophilized
Size	250 µg
Buffer	Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10% Glycerol, pH 7.4
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	<i>Bertholletia excelsa</i> ; Ber e 1
----------	---------------------------------------