



Recombinant Common carp Cyp c 1 (DAG-WT3458)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	12.5 kDa
Format	Lyophilized
Size	250 µg
Buffer	Lyophilized from sterile 20 mM Tris, 200 mM NaCl, pH 7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Common carp; Cyp c 1
----------	----------------------