



# Recombinant *Polistes dominulus* Pol d 5 (DAG-WT3461)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                  |                                                                      |
|------------------|----------------------------------------------------------------------|
| Purity           | > 90% as determined by SDS-PAGE                                      |
| Conjugate        | His                                                                  |
| Applications     | ELISA                                                                |
| Molecular Weight | 36 kDa                                                               |
| Format           | Lyophilized                                                          |
| Size             | 250 µg                                                               |
| Buffer           | Lyophilized from sterile 20 mM Tris, 300 mM NaCl, 10% Glycerol, pH 7 |
| Storage          | Store at -20°C to -80°C                                              |

## BACKGROUND

|          |                                     |
|----------|-------------------------------------|
| Keywords | <i>Polistes dominulus</i> ; Pol d 5 |
|----------|-------------------------------------|