



Recombinant *Oryctolagus cuniculus* Ory c 3 (DAG-WT3464)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	18.3 kDa
Format	Lyophilized
Size	250 µg
Buffer	Lyophilized from sterile 20 mM Tris, 300 mM NaCl, 10% Glycerol, pH 7.5
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	<i>Oryctolagus cuniculus</i> ; Ory c 3
----------	----------------------------------------