



## Recombinant CPP Antigen (DAG-WT3467)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Species	Human
Purity	> 95%
Conjugate	N/A
Applications	ELISA, CLIA, LFIA
Molecular Weight	28kDa
Format	Liquid
Size	100 µg
Buffer	25mM Tris-HCl, 0.02% SDS, pH8.5
Preservative	None
Storage	Store at 2°C to 8°C short term (up to 1 week) or frozen below -20°C

### BACKGROUND

Keywords	Copeptin; CPP; AVP
----------	--------------------