



Recombinant Citrus leprosis virus C p29 (DAG-WT3469)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	Immunogen
Molecular Weight	29.2 kDa
Format	Liquid
Size	1 mg
Buffer	Tris-HCl buffer
Preservative	None
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Citrus leprosis virus C; p29; Citrus sinensis
----------	---

GENE INFORMATION

UniProt ID	Q1KZ58
------------	--------