



Recombinant *Sesamum indicum* Ses I 1 [His] (DAG-WT3470)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	N/A
Molecular Weight	16.4 kDa
Format	Liquid
Size	250 µg
Buffer	PBS buffer, pH 7.4
Storage	Store at -20°C

BACKGROUND

Keywords	Sesamum indicum; Ses I 1
----------	--------------------------