



## Purified Sarcosine oxidase (DAG-WT3471)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	Sarcosine oxidase is obtained from E. coli using proprietary methods. The product is supplied on enzyme activity (KU; The amount of enzyme which produces 1 umol of formaldehyde per min at 37°C and pH 7.7
<b>Conjugate</b>	N/A
<b>Applications</b>	N/A
<b>Bio-activity</b>	>20 U/mg
<b>Format</b>	Lyophilized
<b>Size</b>	10 KU
<b>Storage</b>	Store at -20°C

### BACKGROUND

<b>Keywords</b>	Sarcosine Oxidase; SOX; soxA
-----------------	------------------------------