



# Recombinant Mouse Beta-NGF Protein (DAG-WT3479)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                  |  |
|------------------|--|
| Species          | Mouse  |
| Purity           | > 95% as determined by SDS-PAGE                        |
| Conjugate        | N/A  |
| Applications     | N/A  |
| Molecular Weight | 13.5 kDa   |
| Format           | Lyophilized  |
| Size             | 20 µg, 100 µg, 1 mg                                    |
| Buffer           | Lyophilized from sterile 150mM NaCl, 50mM NaAc, pH 5.5 |
| Storage          | Store at -20°C to -80°C                                |

## BACKGROUND

|          |   |
|----------|---|
| Keywords | Beta-nerve growth factor isoform B precursor; Beta-NGF; NGF |
|----------|---|

## GENE INFORMATION

|            |                |
|------------|----------------|
| UniProt ID | NP_001106168.1 |
|------------|----------------|