



Recombinant Human Beta-NGF Protein (DAG-WT3480)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	N/A
Molecular Weight	24.8 kDa
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from sterile PBS, pH 7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Beta-nerve growth factor isoform B precursor; Beta-NGF; NGF
----------	---

GENE INFORMATION

UniProt ID	NP_002497.2
------------	-------------