



# Recombinant Human Beta-NGF Protein (DAG-WT3481)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Human
Purity	> 85% as determined by SDS-PAGE
Conjugate	N/A
Applications	N/A
Molecular Weight	25 kDa
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from sterile 0.01 % TFA, 40 % acetonitrile
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	Beta-nerve growth factor isoform B precursor; Beta-NGF; NGF
----------	---

## GENE INFORMATION

UniProt ID	NP_002497.2
------------	-------------