



Recombinant Human GLP1R Protein [hFc] (DAG-WT3483)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	hFc
Applications	N/A
Molecular Weight	52.3 kDa
Format	Lyophilized
Size	100 µg, 1 mg
Buffer	Lyophilized from sterile PBS, pH 7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords Glucagon-like peptide 1 receptor; GLP1R; GLP-1 receptor

GENE INFORMATION

UniProt ID NP_002053.3