



# Recombinant Human Tau 412 (1N4R) Protein [His] (DAG-WT3484)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	N/A
Molecular Weight	43.8 kDa
Format	Lyophilized
Size	100 µg, 1 mg
Buffer	Lyophilized from sterile 50mM Tris, 150mM NaCl, 1mM tcep, 1mM EDTA, pH 7.0
Storage	Store at -20°C to -80°C

## BACKGROUND

**Keywords** Microtubule-associated protein tau; MAPT; Neurofibrillary tangle protein; Tau 412; 1N4R

## GENE INFORMATION

**UniProt ID** NP\_001116539.1