



Recombinant Human Tau 441 (2N4R) Protein [His] (DAG-WT3486)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	N/A
Molecular Weight	46.7 kDa
Format	Lyophilized
Size	100 µg, 1 mg
Buffer	Lyophilized from sterile 50mM Tris, 150mM NaCl, 1mM tcep, 1mM EDTA, pH 7.0
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords Microtubule-associated protein tau; MAPT; Neurofibrillary tangle protein; Tau 441; 2N4R

GENE INFORMATION

UniProt ID P10636-8