



Recombinant Human Tau Protein [His] (DAG-WT3487)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 85% as determined by SDS-PAGE
Conjugate	His
Applications	N/A
Molecular Weight	38.7 kDa
Format	Lyophilized
Size	50 µg
Buffer	Lyophilized from sterile PBS, pH 7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Microtubule-associated protein tau isoform 4; Tau
----------	---

GENE INFORMATION

UniProt ID	NP_058525.1
------------	-------------