



Recombinant Human GSK3B Protein [His] (DAG-WT3488)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	N/A
Molecular Weight	50.4 kDa
Format	Liquid
Size	50 µg
Buffer	20mM Tris, 500mM NaCl, pH 7.4, 25% glycerol, 0.5mM PMSF, 0.5mM EDTA
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords Glycogen synthase kinase-3 beta isoform 1; GSK3B

GENE INFORMATION

UniProt ID NP_002084.2