



Recombinant Mouse GSK3B Protein [His] (DAG-WT3489)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Mouse
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	N/A
Molecular Weight	47 kDa
Format	Liquid
Size	50 µg
Buffer	20mM Tris, 500mM NaCl, 25% glycerol, 0.2mM DTT, pH 7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords Glycogen synthase kinase-3 beta isoform 1; GSK3B

GENE INFORMATION

UniProt ID Q9WV60