



# Recombinant Human Beta-APP39 Protein [His, GST] (DAG-WT3491)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                         |   |
|-------------------------|---|
| <b>Species</b>          | Human   |
| <b>Purity</b>           | > 85% as determined by SDS-PAGE                               |
| <b>Conjugate</b>        | His, GST  |
| <b>Applications</b>     | N/A   |
| <b>Molecular Weight</b> | 34 kDa  |
| <b>Format</b>           | Liquid  |
| <b>Size</b>             | 100 µg  |
| <b>Buffer</b>           | 50 mM Tris, 100 mM NaCl, 20% glycerol, 0.05% Tween 20, pH 9.5 |
| <b>Storage</b>          | Store at -20°C to -80°C                                       |

## BACKGROUND

|                 |   |
|-----------------|---|
| <b>Keywords</b> | Amyloid-beta precursor protein; APP; APP-beta; Beta-APP39 |
|-----------------|---|

## GENE INFORMATION

|                   |          |
|-------------------|----------|
| <b>UniProt ID</b> | P05067-1 |
|-------------------|----------|